

**WHAT
FOR
WHAT**

AIR-FLOW® POWDERS

EMS⁺

AIR - FLOW® POWDER		
ORIGINAL AIR - FLOW® THERAPY → SUPRAGINGIVAL	CLINICAL APPLICATIONS	AIR - FLOW® HANDPIECES
 AIR - FLOW® POWDER CLASSIC > Sodium bicarbonate > Average grain size ~65µm > Available in 6 different tastes 300 g bottle	→ REMOVAL OF BIOFILM AND HEAVY DISCOLORATIONS → SURFACE PREPARATION BEFORE DENTAL TREATMENTS	 AIR - FLOW®
 AIR - FLOW® POWDER CLASSIC COMFORT > Sodium bicarbonate > Average grain size ~40µm 250 g bottle	→ REMOVAL OF BIOFILM AND HEAVY DISCOLORATIONS WITH OPTIMIZED PATIENT COMFORT → SURFACE PREPARATION BEFORE DENTAL TREATMENTS	
 AIR - FLOW® POWDER SOFT > Glycine-based > Average grain size ~65µm 200 g bottle	→ REMOVAL OF BIOFILM AND MEDIUM/LIGHT DISCOLORATIONS → TOOTH POLISHING	
→ SUB - AND SUPRAGINGIVAL		
 AIR - FLOW® POWDER PLUS > Erythritol-based > Average grain size ~14µm > Contains chlorhexidine (0,3%) 100 g bottle	→ REMOVAL OF BIOFILM AND MEDIUM/LIGHT DISCOLORATIONS WITH OPTIMIZED PATIENT COMFORT → TOOTH POLISHING → REMOVAL OF BIOFILM FROM POCKETS DOWN TO 4mm ON NATURAL TEETH AND IMPLANTS → REMOVAL OF BIOFILM FROM DEEP PERIODONTAL POCKETS, FROM 4 TO 10mm ON NATURAL TEETH AND IMPLANTS	 AIR - FLOW® AIR - FLOW® PLUS
→ SUBGINGIVAL		
 AIR - FLOW® POWDER PERIO > Glycine-based > Average grain size ~25µm 120 g bottle	→ REMOVAL OF BIOFILM FROM POCKETS DOWN TO 4mm ON NATURAL TEETH AND IMPLANTS → REMOVAL OF BIOFILM FROM DEEP PERIODONTAL POCKETS, FROM 4 TO 10mm ON NATURAL TEETH AND IMPLANTS	 PERIO - FLOW® AIR - FLOW® PERIO



NEW



NEW

**"I FEEL
GOOD"**

EMS-SWISSQUALITY.COM

EMS
ELECTRO MEDICAL SYSTEMS SA
Chemin de la Vuarpilliére 31
CH-1260 Nyon, Switzerland

Tel. +41 22 99 44 700
Fax +41 22 99 44 701
welcome@ems-ch.com
www.ems-dent.com